


Anatomical Breast


Nagasil Smooth	Nagotex Textured	Tissue Expansion Devices Valve:				
		Volume (ml)	Valve	Base (mm)	Height (mm)	Projection (mm)
TA 200	TA-X 200	200	S	100	105	35
TA 250	TA-X 250	250	S	105	115	40
TA 300	TA-X 300	300	S	110	120	45
TA 350	TA-X 350	350	S	110	125	50
TA 400	TA-X 400	400	S	120	135	45
TA 500	TA-X500	500	S	130	140	55
TA 550	TA-X 550	550	S	135	145	55
TA 650	TA-X 650	650	S	145	155	55
TA 750	TA-X 750	750	S	145	155	65



Anatomical
Tissue
Expanders are
specially
designed to
compliment
CoGel®XA
mammary
implants.


Round Breast

Nagasil Smooth	Nagotex Textured	Tissue Expansion Devices Valve:			
		Volume (ml)	Valve	Base (mm)	Projection (mm)
TRD 150	TRD-X 150	150	M	95	40
TRD 200	TRD-X 200	200	M	100	44
TRD 250	TRD-X 250	250	S	110	45




Cylindrical-Straight

Nagasil Smooth	Nagotex Textured	Tissue Expansion Devices Valve: S - Standard dome 26mm M - Mini dome 14mm			
Product Code	Product Code	Volume (ml)	Valve	length (mm)	Diameter (mm)
TCY-S 100	TCY-S-X 100	100	M	110	35
TCY-S 200	TCY-S-X 200	200	M	140	45
TCY-S 300	TCY-S-X 300	300	S	160	50
TCY-S 400	TCY-S-X 400	400	S	175	60




Cylindrical-Curved

Nagasil Smooth	Nagotex Textured	Tissue Expansion Devices Valve: S - Standard dome 26mm M - Mini dome 14mm				
Product Code	Product Code	Volume (ml)	Valve	length (mm)	Curve width (mm)	Diameter (mm)
TCY-C 125	TCY-C-X 125	125	M	120	61	35
TCY-C 250	TCY-C-X 250	250	S	145	77	45
TCY-C 450	TCY-C-X 450	450	S	180	93	55




Croissant

Nagosil Smooth	Nagotex Textured	Tissue Expansion Devices Valve:					
		Volume (ml)	Valve	A (mm)	B (mm)	C (mm)	D (mm)
TCS 125	TCS-X 125	125	M	100	75	70	48
TCS 250	TCS-X 250	250	S	130	93	100	57
TCS 500	TCS-X 500	500	S	170	127	120	78




Calf

Nagosil Smooth	Nagotex Textured	Tissue Expansion Devices Valve:				
		Volume (ml)	Valve	length (mm)	Base (mm)	Projection (mm)
TCF 60	TCF-X 60	60	S	175	40	18
TCF 150	TCF-X 150	150	S	215	40	27
TCF 235	TCF-X 235	235	S	245	55	33



Vaginal

Nagosil Smooth	Nagotex Textured	Product Code	Product Code	Volume (ml)	Valve	length (mm)	Width (mm)
		TEV 80	TEV-X 80	80ml	2-way check	80	40



Accessories

Product Code	Description
TVS	Standard dome valve 26mm with 80mm tube attachment
TVM	Mini dome valve 14mm with 80mm tube attachment